





# ELISA PLATES



Genaxy ELISA plates can be used for wide applications in the field of therapeutics, diagnostics, molecular biology, pharmaceuticals and biomedical studies.

## Unique Features


- » Manufactured with 99.9% virgin polystyrene.
- » Conform to SBS/ANSI this can be used for manual and automated sample processing.
- » Flat bottom and Round bottom formats with or without lid.
- » High binding detachable and non-detachable plate designs are also available.
- » Well volume: 350ul.
- » Protein-binding capacity of High binding plates can reach 550-690ng IgG/cm<sup>2</sup>.
- » All products are sterilized and individually packed.
- » Certified DNase, RNase, Human DNA and Pyrogen free.




# GENAXY ELISA PLATES





Condensation rings on lid  
prevents cross contamination



Excellent well uniformity



Transparent round bottom surface



Transparent flat bottom surface



FOR MORE PRODUCT  
INFORMATION  
SCAN THE QR CODE



## Microtiter Plates

| Cat. Number     | Description                                                | Qty. / Pack |
|-----------------|------------------------------------------------------------|-------------|
| GEN-MTP-96F-S   | Microtitration Plates, Flat Bottom, Sterile                | 100 Plates  |
| GEN-MTP-96U-S   | Microtitration Plates, Round Bottom, Sterile               | 100 Plates  |
| GEN-MTP-L-S     | Lids for Microtiter Plates, Sterile                        | 100 Plates  |
| GEN-MTP-96F-L-S | Microtitration Plates with Lids, Flat Bottom, Sterile      | 100 Plates  |
| GEN-MTP-96U-L-S | Microtitration Plates with Lids, Round Bottom, Sterile     | 100 Plates  |
| GEN-MTP-HB-96F  | 96 Wells Undetachable Elisa Plate, High Binding            | 100 Plates  |
| GEN-MTPS-HB-96F | 96 Wells 8x12 Strips Detachable Elisa Plate, High Binding. | 100 Plates  |

**Genaxy Scientific Pvt. Ltd.**

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058  
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : [info@genaxy.com](mailto:info@genaxy.com)