


GENAXY GEL LOADING TIPS

Genaxy Gel Loading Tips



Gel Loading tips are designed for easy loading of samples on gels. The tips fit a wide variety of single and multichannel pipettors. These tips are different from standard tips, these gel loading tips are designed for most gel electrophoresis tests. The extended tip design is suitable for agarose gel and protein gel loading operation. Soft tip, no damage to gel, allows you to experiment without worries. All tips are certified DNase, RNase and Pyrogen free.

Unique Features

- » Round gel loading tips.
- » Easy and precise loading of samples
- » Fit most popular single and multichannel pipettors, except Rainin LTS and Thermo Clip-Tip.
- » Made of medical PP, Free of RNase, DNase, DNA, and pyrogens.

Cat. Number	Description	Qty. / Pack
GEN-TGL-320R	10µl MicroFlex Gel Loading Tips, Natural, Racked, Flat Tip 0.4mm, Non-Sterile	10 X 96 Tips
GEN-TGL-320R-S	10µl MicroFlex Gel Loading Tips, Natural, Racked, Flat Tip 0.4mm, Sterile	10 X 96 Tips
GEN-TGL-520R	200µl MicroFlex Gel Loading Tips, Natural, Racked, Round Tip 0.57mm, Non-Sterile	10 X 96 Tips
GEN-TGL-520R-S	200µl MicroFlex Gel Loading Tips, Natural, Racked Round Tip 0.57mm, Sterile	10 X 96 Tips

Genaxy Scientific Pvt. Ltd.