




PC CRYO BOX

- The polycarbonate (PC) construction ensures durability.
- Suitable for use at ultra-low temperatures.
- Features a clear, non-reversible lid for better visualization of vials.
- Includes a numeric grid for easy identification of wells/vials.

Cat. No.	Product Description	Pack Size
GEN-CVR-100-PC	100-Well PC Rack for 1.2/2.0ml tubes.	5

Genaxy Scientific Pvt. Ltd .

155-156, C-Block, Community Centre, Plot 32 to 57, Behind Janak Cinema, Janak Puri, New Delhi-110 058
Phones : +91-11-4761 9999 (100 Lines) Fax : +91-11-2553 4163, E-mail : info@genaxy.com